







# Animal Health Laboratory

## Acceptable Serology Samples

<p>Gold-top serum separator tubes (SST)</p>	<p>Contains a clot activator and serum gel separator</p>	
<p>Mottled red/gray-top serum separator tubes (SST)</p>	<p>Contains a clot activator and serum gel separator</p>	
<p>Red-top tube (plastic)</p>	<p>Contains a clot activator but no serum gel separator</p>	
<p>Red-top tube (glass)</p>	<p>Plain glass tube with no clot activator or serum gel separator</p>	
<p>White-top tube (plastic)</p>	<p>Plain plastic tube with no clot activator or serum gel separator</p>	
<p>Clear-top tube (plastic)</p>	<p>Plain plastic tube with no clot activator or serum gel separator</p>	
<p>Eppendorf tubes</p>	<p>Plain plastic tube with no clot activator or serum gel separator</p>	